

# SMART SENSE

BY DIGI<sup>®</sup>

Pharmaceutical

## Protect Medication Safety & Efficacy and Mitigate Risk.

SmartSense is a comprehensive GxP environmental monitoring solution for medication quality and compliance. Eliminate the gap of visibility by implementing a routine condition monitoring solution. The SmartSense platform enables easy remote management of temperature and humidity conditions with real-time alerting (text/phone/email) and reporting. The centralized cellular-based platform helps you maintain medication efficacy while reducing loss with real-time monitoring that surpasses compliance standards.

## A Comprehensive GxP Sensing-as-a-Service Solution



**Mitigate risk with temperature and humidity monitoring**, enabling real-time visibility and alerting.



**Ensure compliance with NIST**  
NIST / ISO 17025 certified and calibrated sensors.



**Leverage real-time data and analytics** for process improvement while identifying underperforming equipment.



**Cellular connectivity removes IT complexity** and allows for real-time visibility during outages.



**Easy to use and install wire-free sensors** with battery life up to 5 years (AA replaceable batteries).



**Versatile placement** in refrigerators, freezers, and warehouses — and integration with existing systems.

## Automate Compliance and Mitigate Risk

"We required a solution to eliminate the gaps in our routine condition monitoring. SmartSense removes this risk and gives us the consistency and visibility we lack. Our teams now focus more on producing and delivering critical pharmaceutical products to our customers."

Senior Director  
Global Pharmaceutical Manufacturer



NIST



## Simple NIST Recalibration Process

With SmartSense, our 3 step NIST recalibration process is quick and easy:

1. A new NIST calibrated probe is sent to you.
2. Simply replace the old NIST probe.
3. The new NIST probe automatically uploads the new calibration certificate to SmartSense.


It's a painless twist on calibrations that takes seconds, resulting in zero interruption to temperature readings. Never again return a sensor to the manufacturer.



SmartSense NIST Probes

## Made for GxP Environments

SmartSense maintains the highest level of compliance with industry standards to ensure our products and services meet the regulatory requirements critical to our clients. Our solution includes a GxP validated environment.



For more information about SmartSense visit [www.smartsense.co](http://www.smartsense.co) or call +1 (866) 806-2653

**SMART SENSE**  
BY DIGI